

Guidelines for recommending Bambach seat

STEP 1

Select a Seat with or without a Back

Seat with No Back
Seat with a Back

for tasks when seated for short periods eg dentistry, hospitals, medical Imaging
for tasks that require sitting for long periods, eg office workers.

STEP 2

Select Size of seat

Small

Recommended for Women 5ft 4' and below (up to size 10) and children up to 13ys.

Medium

Children 14yrs +,
Male adults up to 5'10" – slim build
Female adults up to 5'10" – standard build

Narrow

Females 5'2" – 5'6" – slim build

The NARROW Seat is essentially the same size seat base as the Medium, however the sides are 'cut away' which allows the user to sit with their legs closer together. Often preferred by Ladies.

Large

Males over 5' 7" (Slender gentlemen may still prefer the Medium) or Larger Females

The LARGE Seat is a larger version of the STANDARD seat, with extra padding , a greater width and length allows for more comfort in the genital area (for males), a more supporting seat for larger individuals.



STEP 3

Select Gas Stem – 4 options

Based on the sitters height, use the guidelines below to determine the correct Gas Stem to be fitted:

GAS STEMS	Customer Height	GAS STEMS	Customer Height
LOW 100mm	4'2"-4'8"	STD 200mm	5'2"-6'2"
MED 140mm	4'8"-5'3"	DRFT 270mm	6'0"-6'5"

STEP 4

Select Castor/Wheel type

Based on the 'type of flooring' the Seat will be used on, use the guidelines below to determine the correct Castor Type:

STANDARD – carpet	REV. FRICTION – lock without weight
SOFT – hard surfaces, vinyl, wood, slate	CHROME HOODED – carpet use
LOCKING – clip locks on Wheel	GLIDES – flat, do not move at all
FRICTION –lock with weight	LARGE – suitable for hospitals


STEP 5

Seat with or without armrest. Please advise if armrest should be positioned **LEFT** or **RIGHT**.

STEP 6

Select colour

Standard Vinyl Range



Final recommendation

By changing to a Bambach Saddle Seat, it is recommended you take time to get used to this new way of being seated. Limit the time you are seated to 20 - 30 minute intervals in the first few days, in order to avoid becoming 'saddle sore'. Please review the following website www.bambach.com.au and the informational DVD regarding 'adjusting your Saddle Seat'.

For more information contact Ampac Dental on 1800 061 343.