

Dental Chairs



Patient Chair and Hydro Group	Friend Up Australia	Friend Up
Continental (Swing arm) OR Hanging Hoses	✓	✓
Lifting capacity	200 Kg	200 Kg
Automatic programs	4	4
Soft Seamless ergonomic upholstery	✓	✓
Double articulated headrest	✓	✓
Faro Operator Light	Maya LED	Edi Non LED
Single arm rest	✓	✓
Foldable right arm rest	✓	Optional
Removable and autoclavable ceramic cuspidor for easy cleaning. Cuspidor rotates 180° for patient easy access.	✓	✓
Automatic Cup fill and bowl rinse	✓	✓
Multi-functional foot control for easy access to chair positions, chair adjustments, exit position, rinse position and handpiece control.	✓	Optional
Foot control for easy access of instrument functions	✓	✓
Emergency cut off switch	✓	✓
Semi Dry Suction compatability	✓	✓
Easily accessible suction filter	✓	✓
Durr place selector valve	✓	✓
Cattani spittoon valve	✓	✓
Clean Water Bottle System	✓	Optional
Nurses Console		
Assistants console complete with Luzzani Triplex syringe, Cattani High and Low volume suction hoses.	✓	✓
Operation of the cup fill and bowl flush (automatic or manual)	✓	✓
Operation of patient chair and instrument functions	✓	Optional
Doctors Instrument Console		
Control panel allows 4 separate chair programs	✓	✓
Instrument positions include Luzzani Triplex syringe, 1 x fibre optic high speed tubing, 1 Bien Air Micromotor Fibre Optic MC3LK (value \$1800+gst) and space for scaler	✓	✓
Second high speed fibre optic line	✓	Optional
Doctors Stools	Optional	Optional
Assistants' Stools with foot ring	Optional	Optional
Colours Choose from 23 contemporary colours	✓	✓
Warranty		
2 years on parts (excluding globes, handpiece tubing, hoses and suction valves) and 1 year labour (Labour at customers' premises only excluding travelling expenses)	✓	✓