




ampacdental
DENTAL EQUIPMENT SPECIALISTS

Nardi dental compressors are manufactured in Italy to deliver high quality, **OIL FREE, CLEAN, DRY AIR.**



25L 73 dB



25L hood 67dB



30L 64 dB; 50L 62 dB



- Oil Free, Dry Air
- Maintenance Free
- 4 years parts warranty
- Made in Italy
- Certified by Workcover NSW
- Suitable for HOT climates



50L 67dB



50L 73 dB



ampac dental

DENTAL EQUIPMENT SPECIALISTS



50L 73 dB



x 2



90L 73 dB



x 3



90L 73 dB CAD\CAM

Oil residue can damage handpieces and fine pneumatic controls in the dental units.

An oil free compressor is vital in protecting the health of your patients, staff and ultimately increase the longevity of your equipment.

tank size	dryer	cylinders	air intake l/min	output l/min @5bar	pressure Bar	power HP	decibels dB	weight Kg	dimension L x W x H cm
25 l	yes	2	350	156	8	2.00	73	44	59 x 43 x 66
25 l hood	yes	2	350	156	8	2.00	67	48	70 x 43 x 78
30 l cabinet	yes	2	350	156	8	2.00	64	79	71 x 44 x 76
50 l W	yes	2	350	156	8	2.00	73	51	84 x 44 x 68
50 l T	yes	4	480	240	8	2.50	73	60	64 x 45 x 89
50 l hood	yes	2	350	156	8	2.00	67	55	70 x 43 x 98
50 l cabinet	yes	2	350	156	8	2.00	62	96	85 x 50 x 82
90 l	yes	4	480	240	8	2.50	73	86	110 x 59 x 81
90 l cad/cam	2	4	480	240*	10	2.50	73	86	110 x 59 x 81

*output @ 8 bar=190 l/min

MADE IN ITALY



ampacdental.com.au • info@ampac-dental.com.au