

With Swinging Instruments Arm:

It's a new generation unit, a combination of the high Swiss technology with Italian production and design. Ergonomic handling, smooth operating and a high quality at very competitive price.



The Hydro Unit :

- Removable and auto-cleavable ceramic cuspidor bowl
- Automatic timer of cup filler and bowl flush.
- Tube of 20" for a direct discharge
- Filters for aspiration and cuspidor bowl are inspected from the outside part.

The Electromechanical Patient Chair "AG03SWF" with motor lifting to low tension

- The headrest is provided with a double articulation
- 4=3 different chair programs positions and plus one additional rising and return to the last work position
- All positions are stored and controlled from the dentist side and part from the assistant's side
- The chair is equipped with 2 safety switches
- No cables or supply hoses on floor

Spécifications Techniques/Détails :

Unit de base Friend easy Art .Nr. SW77520 with delivery arm system:

(Instrumentation From left to right):

Pos.1=1 3-way f syringe,

Pos.2=(free),

Pos.3=(free),

Pos.4=1 Air motor Vortige with integrated spray without light

Pos.5=1 Turbine tubing without light.

- Predisposition of light (Fibre optic) optional
- Extensions options : 2nd Turbine or 2nd air Motor, air scalar AIRSON Or curing light
- Control panel for patient's chair and instruments
- Swivel tray support with one tray



Hydro Group unit :

- Removable and Auto-cleavable ceramic cuspidor bowl
- Can rotate 180° (space saving)
- Cuspidor's drainage valve included
- Tube of 20" for a direct discharge
- Predisposition for Aqua distilled bottle (optional)



Unité Assistante: (Instrumentation de gauche à droite)

- Extensions options to 3 instruments

Pos.1=1 Surgical section,

Pos.2=1 saliva ejector,

Pos.3=(free =optional : 1 optional 3-way-syringe for assistant)



Fauteuil patient AG03SWF:

- 4=3 programs positions, plus additional rising position and return to the last position.
- Left arm rest (Right one optional).
- The Head rest with a double articulation (Bios with mechanical or pneumatic inclination optional).



The pneumatic foot control:

- The foot control for the command of instruments functions only.
- Activation of the instruments with a progressive regulation of micro motors power.
- The spray and chip blower function are included.
- Independent foot pedal for chair movement optional



The operating Lamp: (Unit Mounted)

- Equipped with EDI Light (Alya Led optional)
- Reliable and proved
- 24.000 Lux, 5000 °k



Colors :

- For exact definition colors please refer to the separate Chair's upholstery color table
- Unit base, arm systems, and Hydro unit standard color :(color choose optional)